

Electrosurgery

BipoJET®	Bipolar	Scissors
.....	1168	Non-Irrigating Bipolar
Forceps1169 - 1190	Irrigating Bipolar
Forceps1191 - 1194	Disposable
Bipolar Forceps1195 - 1197	Bipolar
Connecting Cords1198 - 1200	

Electrosurgery

BipoJET® Bipolar Scissors

The versatility of BipoJET Dissecting Bipolar Scissors enables:

Parallel Cutting and Coagulating



Precoagulation



Pinpoint Coagulation



Surface Coagulation



For more information about Bipolar Scissors see our brochure number **DOC456**.
See page 1199-1200 for appropriate connecting cords.



BC691R-BC693R



BC691R-BC693R 1/1

METZENBAUM Bipolar Dissecting Scissors

	Catalog Number	Length (in.)
BC691	X-	4 3/4
R	FINE	5 1/2
BC692	X-	6 1/4
R	FINE	
BC693	X-	
R	FINE	



BC694R-BC696R



BC694R-BC696R 1/1

STEVENS Bipolar Dissecting Scissors

Catalog Number	Tip	Length (in.)
BC694	TENOTOM	4 3/4
R	Y	5 1/2
BC695	TENOTOM	6 1/4
R	Y	
BC696	TENOTOM	
R	Y	



BC685R-BC690R

1/1



BC685R 1/1




BC686R 1/1

METZENBAUM Bipolar Dissecting Scissors

	Catalog Number	Length (in.)
BC685	STANDARD	7
R	STANDARD	9
BC687	STANDARD	11
R	FINE	7
BC689	FINE	9
R	FINE	11

BC686
R
BC688
R
BC690
R



YASARGIL Bipolar Forceps Features:

- High precision, tissue preserving preparation instrument
- Stronger tension of the forceps allows for blunt spreading

of the tissue.

Tip Size Reference Chart

0.2 mm		0.8 mm	
0.25 mm		0.9 mm	
0.4 mm		1.0 mm	
0.5 mm		1.3 mm	
0.6 mm		1.5 mm	
0.7 mm		2.0 mm	

20 mm, 3/4"



35 mm, 1 3/8"



55 mm, 2 1/4"



1/1
GK746R, GK760R, GK764R

1/1
GK747R, GK761R, GK765R

1/1
GK748R, GK762R, GK766R

1/1
GK749R, GK763R, GK767R

YASARGIL Style Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK746R	0.4	BAYONET	INS*	NON-IRRIG*	3/4	5 1/2
GK747R	0.7	BAYONET	INS	NON-IRRIG	3/4	5 1/2
GK748R	1.0	BAYONET	INS	NON-IRRIG	3/4	5 1/2
GK749R	1.3	BAYONET	INS	NON-IRRIG	3/4	5 1/2
GK760R	0.4	BAYONET	INS	NON-IRRIG	1 3/8	6 1/8
GK761R	0.7	BAYONET	INS	NON-IRRIG	1 3/8	6 1/8
GK762R	1.0	BAYONET	INS	IRRIG	1 3/8	6 1/8
GK763R	1.3	BAYONET	INS	IRRIG	1 3/8	6 1/8
GK764R	0.4	BAYONET	INS	IRRIG	2 1/4	6 7/8
GK765R	0.7	BAYONET	INS	NON-IRRIG	2 1/4	6 7/8
GK766R	1.0	BAYONET	INS	NON-IRRIG	2 1/4	6 7/8
GK767R	1.3	BAYONET	INS	NON-IRRIG	2 1/4	6 7/8

* Key:
INS = Insulating
NON-IRRIG = Non-Irrigating
IRRIG = Irrigating

NON-IRRIG


Electrosurgery

Non-Irrigating Bipolar Forceps

YASARGIL Style Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK768R	0.4	BAYONET	INS	NON-IRRIG	3	7 3/4
GK769R	0.7	BAYONET	INS	NON-IRRIG	3	7 3/4
GK770R	1.0	BAYONET	INS	NON-IRRIG	3	7 3/4
GK771R	1.3	BAYONET	INS	NON-IRRIG	3	7 3/4
GK772R	0.4	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK773R	0.7	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK774R	1.0	BAYONET	INS	IRRIG	3 3/4	8 1/2
GK775R	1.3	BAYONET	INS	IRRIG	3 3/4	8 1/2

NON-IRRIG



1/1
GK768R, GK772R

1/1
GK769R, GK773R

1/1
GK770R, GK774R



1/1
GK771R, GK775R

ELECTRO SURGERY

MILLIMETERS

YASARGIL Style Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK788R	0.4	BAYONET	INS	NON-IRRIG	4 1/2	9 1/4
GK789R	0.7	BAYONET	INS	NON-IRRIG	4 1/2	9 1/4
GK790R	1.0	BAYONET	INS	NON-IRRIG	4 1/2	9 1/4
GK791R	1.3	BAYONET	INS	NON-IRRIG	4 1/2	9 1/4
GK793R	0.7	BAYONET	INS	NON-IRRIG	5 1/4	10



-  1/1
GK788R
-  1/1
GK789R, GK793R
-  1/1
GK790R
-  1/1
GK791R







Electrosurgery

Non-Irrigating Bipolar Forceps

YASARGIL Style Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK777R	45° UP	0.7	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK780R	90° UP	0.4	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK781R	90° UP	0.7	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2
GK785R	90° DOWN	0.7	BAYONET	INS	NON-IRRIG	3 3/4	8 1/2



YASARGIL Style Bipolar Forceps


Catalog Number	Tip	Tip Size (mm)	Working Length (in.)	Length (in.)
GK831R	RING	2/	4 1/2	9 1/4
GK832R	RING	3	4 1/2	9 1/4
GK833R	RING	3/	4 1/2	9 1/4

4
4/
5

Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK660R	0.6	BAYONET	INS	NON-IRRIG	1 1/2	6 1/4
GK670R	0.6	BAYONET	INS	NON-IRRIG	3 1/8	7 7/8
GK671R	0.9	BAYONET	INS	NON-IRRIG	3 1/8	7 7/8
GK675R	0.6	BAYONET	INS	NON-IRRIG	4	8 3/4
GK677R	0.6	BAYONET	INS	NON-IRRIG	4 3/4	9 1/2
GK679R	0.9	BAYONET	INS	NON-IRRIG	4 3/4	9 1/2

IRRIG



Electrosurgery

Non-Irrigating Bipolar Forceps

Bayonet Forceps, Rose Gold, Insulated

■ Rose Gold plated tips reduce tissue adhesion



Non-Stick Rose Gold forceps

The Rose Gold forceps enable heat to be quickly discharged away from the tip.



Standard forceps

The heat cannot be discharged fast enough. This leads to overheating and sticking of tissue.

Rose Gold Bipolar Forceps, Bayonet-Shaped

Catalog Number	Tip Size (mm)	Style	Insulation	Irrigation	Working Length (in.)	Length (in.)
US403	0.	HARDY	IN	NON-IRRG	2	6
US404	5	HARDY	S	NON-IRRG	2	6
US405	0.	HARDY	IN	NON-IRRG	2	6
US406	7	HARDY	S	NON-IRRG	2	6
US407	1.	HARDY	IN	NON-IRRG	2	6
US605	0	SCOVILLE	S	NON-IRRG	2 3/4	7 3/4
US606	1.	SCOVILLE	IN	NON-IRRG	2 3/4	7 3/4
US607	5	SCOVILLE	S	NON-IRRG	2 3/4	7 3/4
US608	2.	SCOVILLE	IN	NON-IRRG	2 3/4	7 3/4
US609	0	SCOVILLE	S	NON-IRRG	2 3/4	7 3/4
US600	0.	CUSHING	IN	NON-IRRG	3 1/2	7 1/2
US601	5	CUSHING	S	NON-IRRG	3 1/2	7 1/2
US602	0.	CUSHING	IN	NON-IRRG	3 1/2	7 1/2
US603	7	CUSHING	S	NON-IRRG	3 1/2	7 1/2
US604	1.	CUSHING	IN	NON-IRRG	3 1/2	7 1/2
US610	0	HARDY	S	NON-IRRG	4 1/2	8 1/2
US611	1.	HARDY	IN	NON-IRRG	4 1/2	8 1/2
US612	5	HARDY	S	NON-IRRG	4 1/2	8 1/2
US613	2.	HARDY	IN	NON-IRRG	4 1/2	8 1/2
US614	0	HARDY	S	NON-IRRG	4 1/2	8 1/2
US615	0.	HARDY	IN	NON-IRRG	5	9
US616	5	HARDY	S	NON-IRRG	5	9
US617	0.	HARDY	IN	NON-IRRG	5	9
US618	7	HARDY	S	NON-IRRG	5	9
US619	1.	HARDY	IN	NON-IRRG	5	9
US530	0	HARDY	S	NON-IRRG	6	10
US620	1.	HARDY	IN	NON-IRRG	6	10
US621	5	HARDY	S	NON-IRRG	6	10
US173	2.	HARDY	IN	NON-IRRG	6	10
US531	0	HARDY	S	NON-IRRG	6	10



US173-US621

US403, US530, US600, US610, US615

US404, US601, US611, US616, US621

US405, US602, US612, US617

US173, US406, US603, US613, US618

US407, US531, US604, US614, US619

0. IN
5 S
0. IN
7 S
1. IN
0 S
1. IN
5 S
2. IN
0 S
0. IN
5 S
0. IN
7 S
1. IN
0 S
1. IN
5 S

ELECTRO SURGERY

Rose Gold Microform Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US526	0.5	BAYONET	INS	NON-IRRG	1 3/8	6
US527	0.7	BAYONET	INS	NON-IRRG	1 3/8	6
US528	1.0	BAYONET	INS	NON-IRRG	1 3/8	6
US529	1.3	BAYONET	INS	NON-IRRG	1 3/8	6
US570	0.5	BAYONET	INS	NON-IRRG	2 3/8	7
US571	0.7	BAYONET	INS	IRRG NON-	2 3/8	7
US572	1.0	BAYONET	INS	IRRG NON-	2 3/8	7
US573	1.3	BAYONET	INS	IRRG NON-	2 3/8	7
US574	0.5	BAYONET	INS	IRRG	3	7 3/4
US575	0.7	BAYONET	INS	NON-IRRG	3	7 3/4
US576	1.0	BAYONET	INS	NON-IRRG	3	7 3/4
US577	1.3	BAYONET	INS	NON-IRRG	3	7 3/4
US578	0.5	BAYONET	INS	NON-IRRG	3 3/4	8 1/2
US579	0.7	BAYONET	INS	NON-IRRG	3 3/4	8 1/2
US580	1.0	BAYONET	INS	IRRG NON-	3 3/4	8 1/2
US581	1.3	BAYONET	INS	IRRG NON-	3 3/4	8 1/2
US582	0.5	BAYONET	INS	IRRG	4 1/2	9 1/4
US583	0.7	BAYONET	INS	IRRG	4 1/2	9 1/4
US584	1.0	BAYONET	INS	NON-IRRG	4 1/2	9 1/4
US585	1.3	BAYONET	INS	NON-IRRG	4 1/2	9 1/4
US586	0.5	BAYONET	INS	NON-IRRG	5 1/2	10 1/4
US587	0.7	BAYONET	INS	NON-IRRG	5 1/2	10 1/4
US588	1.0	BAYONET	INS	IRRG NON-	5 1/2	10 1/4
US589	1.3	BAYONET	INS	IRRG NON-	5 1/2	10 1/4



IRRG
NON-
IRRG
NON-
IRRG



ELECTRO SURGERY

Electrosurgery

Non-Irrigating Bipolar Forceps

Rose Gold Reverse Action Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US547	0.7	BAYONET	INS	NON-IRRG	6	10
US548	1.0	BAYONET	INS	NON-IRRG	6	10



Rose Gold Minimally Invasive (MIS) Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US541	STR 0.7	BAYONET		INS	NON-IRRG	5 1/4	9 1/2
US542	STR 1.2	BAYONET		INS	NON-IRRG	5 1/4	9 1/2
US543	STR 0.7	BAYONET		INS	NON-IRRG	6 1/4	10 1/2
US544	STR 1.2	BAYONET		INS	NON-IRRG	6 1/4	10 1/2
US545	45° UP0.7	BAYONET		INS	NON-IRRG	5 1/4	9 1/2
US546	45° UP0.7	BAYONET		INS	IRRG	6 1/4	10 1/2




Rose Gold SEMKIN Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US413	1.5	STR	INS	NON-IRRG	6



SCOVILLE Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US534	1.0	BAYONET	INS	NON-IRRIG	2 3/4	7 3/4
US301	1.5	BAYONET	INS	NON-IRRIG	2 3/4	7 3/4
US300	1.5	BAYONET	NON-INS	NON-IRRIG	2 3/4	7 3/4



GERALD Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US342	1.0	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2
US307	1.0	BAYONET	NON-INS	NON-IRRIG	3 1/2	7 1/2



CUSHING Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US308	0.7	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2
US182	1.0	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2
US183	1.5	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2
US305	2.0	BAYONET	INS	NON-IRRIG	3 1/2	7 1/2



Electrosurgery


Non-Irrigating Bipolar Forceps

HARDY Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US310	0.5	BAYONET	INS	NON-IRRG	4 1/2	8 1/2
US536	0.7	BAYONET	INS	NON-IRRG	4 1/2	8 1/2
US184	1.0	BAYONET	INS	NON-IRRG	4 1/2	8 1/2
US185	1.5	BAYONET	INS	NON-IRRG	4 1/2	8 1/2
US186	2.0	BAYONET	INS	NON-IRRG	4 1/2	8 1/2
US187	0.5	BAYONET	INS	NON-IRRG	6	10
US188	1.0	BAYONET	INS	IRRG	6	10
US189	1.5	BAYONET <td INS	IRRG	6	10	




US184-US189, US310, US536




US187, US310



US536



US184, US188



US185, US189



US186

HARDY Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US309	0.5	BAYONET	NON-INS	NON-IRRG	4 1/2	8 1/2



US309



US309

ELECTRO SURGERY

MILLIMETERS

BARRÉ BipoJET® Nerve/Tissue Sparing Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
GK685R	3.0	BAYONET	INS	NON-IRRIG	10 1/4



Strong Tension Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US314	0.5	BAYONET	NON-INS	NON-IRRIG	2 3/8	6 1/4
US315	0.7	BAYONET	NON-INS	NON-IRRIG	2 3/8	6 1/4
US316	1.0	BAYONET	NON-INS	NON-IRRIG	2 3/8	6 1/4
US317	1.3	BAYONET	NON-INS	NON-IRRIG	2 3/8	6 1/4
US318	0.5	BAYONET	NON-INS	NON-IRRIG	3 1/4	7 1/8
US319	0.7	BAYONET	NON-INS	NON-IRRIG	3 1/4	7 1/8
US320	1.0	BAYONET	NON-INS	NON-IRRIG	3 1/4	7 1/8
US321	1.3	BAYONET	NON-INS	NON-IRRIG	3 1/4	7 1/8
US322	0.5	BAYONET	INS	NON-IRRIG	4	7 7/8
US323	0.7	BAYONET	NON-INS	IRRIG	4	7 7/8
US324	1.0	BAYONET	NON-INS	IRRIG	4	7 7/8
US325	1.3	BAYONET	NON-INS	IRRIG	4	7 7/8




NON-INS
NON-IRRIG



Electrosurgery

Non-Irrigating Bipolar Forceps



CASPAR Style Bipolar Forceps:

- Developed for precision neurosurgery with a special focus on spine surgery.
- Parallel guiding handle for exact tip alignment.

CASPAR Style Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK899	0.	BAYONE	IN	NON-IRRG	2 3/4	6 1/2
R	2	T	S	NON-IRRG	2 3/4	6 1/2
GK900	0.	BAYONE	IN	NON-IRRG	2 3/4	6 1/2
R	5	T	S	NON-IRRG	2 3/4	6 1/2
GK910	1.	BAYONE	IN	NON-IRRG	4	7 3/4
R	0	T	S	NON-IRRG	4	7 3/4
GK920	2.	BAYONE	IN	NON-IRRG	4	7 3/4
R	0	T	S	NON-IRRG	4	7 3/4
GK929	0.	BAYONE	IN	NON-IRRG	5	8 3/4
R	2	T	S	NON-IRRG	5	8 3/4
GK930	0.	BAYONE	IN	IRRG	5	8 3/4
R	5	T	S	NON-IRRG	5	8 3/4




GK940	1.	BAYONE	IN	IRRG		
R	0	T	S			
GK950	2.	BAYONE	IN			
R	0	T	S			
GK959	0.	BAYONE	IN			
R	2 1/2	T	S			
GK960	0, 1	BAYONE	IN			
R	5	T	S			
GK920R, GK950R, GK9580R						

GK970 1. BAYONE IN

CASPAR Style Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Shape	Tip Size (mm)	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK976R	90° UP	BAYONET	1.0	INS	NON-IRRG	4 3/4	8 1/2
GK978R	90° DOWN	BAYONET	1.0	INS	NON-IRRG	4 3/4	8 1/2
GK972R	45° UP	BAYONET	1.0	INS	NON-IRRG	5	8 3/4
GK974R	45° DOWN	BAYONET	1.0	INS	NON-IRRG	5	8 3/4




Electrosurgery

Non-Irrigating Bipolar Forceps

SINTRAM Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK644RSTR		0.7 x 6	BAYONET	INS	NON-IRRIG	2 3/8	6 1/4
GK645RSTR		1.0 x 6	BAYONET	INS	NON-IRRIG	3 1/4	7 1/4
GK646RSTR		1.0 x 6	BAYONET	INS	NON-IRRIG	4	7 7/8
GK614R	STR	0.7 x 23	BAYONET	INS	NON-IRRIG	2 3/8	6 1/4
GK615R	STR	1.0 x 23	BAYONET	INS	NON-IRRIG	3 1/4	7 1/4
GK616R	STR	1.0 x 23	BAYONET	INS	NON-IRRIG	4	7 7/8
GK604R	STR	0.9 x 23	STR	INS	IRRIG NON-	-	6 1/4
GK605R	STR	0.9 x 23	STR	INS	IRRIG NON-	-	7 1/4
GK606R	STR	0.9 x 23	STR	INS	IRRIG NON-	-	7 7/8
GK624R	ANG 35°	0.9 x 23	STR	INS	NON-IRRIG	-	6 1/4
GK625R	ANG 35°	0.9 x 23	STR	INS	NON-IRRIG	-	7 7/8




NON-IRRIG

ELECTRO SURGERY

MILLIMETERS

MALIS Style Titanium Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US331	STR	0.	BAYONE	NON-INS	NON-IRRG	3 3/4	8 3/4
US332	STR	5	T	NON-INS	NON-IRRG	3 3/4	8 3/4
US330	STR	1.	BAYONE	NON-INS	NON-IRRG	3 3/4	8 3/4
US333	ANG-UP	0	T	NON-INS	NON-IRRG	3 3/4	8 3/4
US175	DOWN-ANG	0.5	BAYONE	NON-INS	NON-IRRG	3 3/4	8 3/4



MALIS Style Rose Gold Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US161	STR	0.5	BAYONET	INS	NON-IRRG	3 3/4	8 3/4
US162	STR	1.0	BAYONET	INS	NON-IRRG	3 3/4	8 3/4
US163	STR	1.5	BAYONET	INS	NON-IRRG	3 3/4	8 3/4
US164	ANG-UP	0.5	BAYONET	INS	NON-IRRG	3 3/4	8 3/4
US165	DOWN-ANG	0.5	BAYONET	INS	NON-IRRG	3 3/4	8 3/4



Electrosurgery

Non-Irrigating Bipolar Forceps

Titanium Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation NON-INS	Irrigation NON-IRR	Length (in.)
US32	STR	0.5	ST	INS	IRR	6
6	STR	0.5	R	INS	IRR	7 1/8
US32	CVD	0.5	ST	NON-INS	NON-IRR	6
8	CVD	0.5	R	NON-INS	NON-IRR	7 1/8



US326, US328



US327, US329

US32			ST			
7			R			
US32			ST			
			R			

Micro Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Irrigation NON-INS	Insulation NON-IRR	Length (in.)
US34	STR	0.7	ST	INS	IRR	6
5	STR	0.7	R	INS	IRR	7 1/8
US34	CVD	0.7	ST	NON-INS	NON-IRR	6
7	CVD	0.7	R	NON-INS	NON-IRR	7 1/8



US345, US347



US346, US348

US34			ST			
6			R			
US34			ST			
8			R			

Classic Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation NON-IRR	Length (in.)
GK700	STR	0.6	ST	IN	IRR	6 1/4
R	STR	0.9	R	S	IRR	6 1/4
GK701	STR	0.9	ST	IN	IRR	7 7/8
R	ANG-DOWN	0.9	R	S	IRR	6 1/4
GK703	ANG-DOWN	0.9	ST	IN	IRR	7 7/8



GK700R-GK703R



GK700R, GK701R, GK703R



GK711R-GK713R



GK711R, GK713R

R	DOWN		R	S		
GK711			ST	IN		
R			R	S		

ELECTRO SURGERY

MILLIMETERS




Classic Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
GK690R	STR	0.	ST	IN	NON-IRRG	4 3/4
GK720R	STR	4	R	S	NON-IRRG	6 1/4
GK722R	STR	0.	ST	IN	NON-IRRG	6 1/4
GK730R	STR	6	R	S	NON-IRRG	7 7/8
GK723R	ANG-DOWN	0.	ST	IN	NON-IRRG	6 1/4

9 R S
0. ST IN
9 R S

0. ST IN



Jewelers Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US337	STR	0.4	STR	INS	NON-IRRG	4
US335	STR	0.4	STR	NON-INS	NON-IRRG	4
US159	CVD	0.4	STR	INS	NON-IRRG	4
US336	CVD	0.4	STR	NON-INS	NON-IRRG	4




TENZEL Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US178	CVD	0.4	STR	INS	NON-IRRG	4
US338	CVD	0.4	STR	NON-INS	NON-IRRG	4



Electrosurgery

Non-Irrigating Bipolar Forceps



US339-US340

ADSON Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US340	STR	1.0	STR	INS	NON-IRRG	5
US339	STR	1.0	STR	NON-INS	NON-IRRG	5



US179-US344

SEMKIN Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US179	0.7	STR	INS	NON-IRRG	5 1/2
US311	0.7	STR	NON-INS	NON-IRRG	5 1/2
US344	2.0	STR	INS	NON-IRRG	6
US196	2.0	STR	NON-INS	NON-IRRG	6



US197-US304

CUSHING Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US197	0.7	STR	INS	NON-IRRG	6
US302	0.7	STR	NON-INS	NON-IRRG	6
US304	1.5	STR	INS	NON-IRRG	7
US303	1.5	STR	NON-INS	NON-IRRG	7

GERALD Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Shape	Insulation	Irrigation	Length (in.)
US341	STR	STR	INS	NON-IRRIG	7
US306	STR	STR	NON-INS	NON-IRRIG	7



US306-US341

Rose Gold Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Style	Insulation	Irrigation	Length (in.)
US171	1.0	STR	ADSON	INS	NON-IRRIG	5
US172	0.4	STR	JEWELERS	INS	NON-IRRIG	4



US171



US171



US172



US172

MALIS Style Titanium Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US176	0.25	STR	NON-INS	NON-IRRIG	6
US177	0.5	STR	NON-INS	NON-IRRIG	6



US176-US177



US176-US177

Micro and Laboratory Forceps, Stainless Steel

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
GK596R	STR	0.3	STR	INS	NON-IRRIG	5 3/4
GK595R	STR	0.6	STR	INS	NON-IRRIG	4 3/4
GK597R	CVD	0.6	STR	INS	NON-IRRIG	4 1/4



GK595R-GK596R



GK595R



GK596R



GK597R



GK597R



Electrosurgery

Non-Irrigating Bipolar Forceps

BiProtect®

Ceramic Insulated Tip, Insulated

■ Ceramic tip insulation helps ensure that heat concentrates on the inner and distal part of the forceps

- Less collateral damage and reduced thermal spread
- Partially insulated tip allows for precise and safe surface coagulation

YASARGIL BiProtect Ceramic Insulated Tip Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GL470R	STR	0.4	BAYONET	INS	NON-IRRG	3 1/2	7 3/4
GL471R	STR	0.8	BAYONET	INS	NON-IRRG	3 1/2	7 3/4
GL480R	STR	0.4	BAYONET	INS	NON-IRRG	4 1/2	8 1/2
GL481R	STR	0.8	BAYONET	INS	NON-IRRG	4 1/2	8 1/2



GL470R-GL481R



GL470R, GL480R



GL471R, GL481R

Transnasal Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK719R	1.0	ANG	INS	NON-IRRG	5 1/8	8 1/2



GK719R



GK719R

130 mm, 5 Pcs

MIN Instruments, Pivot-Point Bipolar Forceps, Transnasal

Catalog Number	Tip	Tip Size (mm)	Shape	Working Length (in.)	Length (in.)
GK822	STR	0.	BAYONE	3 3/4	8 1/2
R	STR	7	T	3 3/4	8 1/2
GK824	STR	1.	BAYONE	5 1/4	10
R	STR	0	T	5 1/4	10
GK826	CVD	0.	BAYONE	3 3/4	8 1/2
R	CVD	7	T	3 3/4	8 1/2
GK828	CVD	1.	BAYONE	5 1/4	10
R	CVD	0	T	5 1/4	10



Special pin between the branches opens the tip of the forceps by additional compression of the handle, allowing secure coagulation in narrow and deep seated surgical fields.

95 mm, 3 7/8"

GK822R-GK824R

135 mm, 5 3/8"


GK826R-GK828R


95 mm, 3 7/8"

GK823R-GK825R

135 mm, 5 3/8"

GK827R-GK829R

 1/1
GK822R, GK823R, GK826R, GK827R

 1/1
GK824R, GK825R, GK828R, GK829R

GK823 0. BAYONE
R 7 T
GK825 1. BAYONE
0 T
GK827 0. BAYONE
R 7 T
GK829 1. BAYONE
R 0 T





Electrosurgery

Non-Irrigating Bipolar Forceps

LANDOLT Bipolar Forceps, Stainless Steel

Catalog Number	Tip	Insulation	Irrigation	Description	Working Length (in.)
GK560R	ANG 90°	INS	NON-IRRIG	COAGULATION FORCEPS F/ HYPOPHYSECTOMY	6
GK561R	90°	-	-	SPARE INNER PART F/ GK560R	
GK580R	ANG 120°	INS	NON-IRRIG	COAGULATION FORCEPS F/ HYPOPHYSECTOMY	6
GK581R	120°	-	-	SPARE INNER PART F/ GK580R	



ELECTRO SURGERY

MILLIMETERS

■ Irrigating coagulation forceps facilitate gentle coagulation without tissue adhesion.

Rose Gold Irrigating Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US622	0.	BAYONE	IN	IRR	3 1/2	7 1/2
US623	5	T	S	G	3 1/2	7 1/2
US624	0.	BAYONE	IN	IRR	3 1/2	7 1/2
US625	7	T	S	G	3 1/2	7 1/2
US626	1.	BAYONE	IN	IRR	3 1/2	7 1/2
US627	0	T	S	G	3	7 3/4
US628	1.	BAYONE	IN	IRR	3	7 3/4
US629	5	T	S	G	3	7 3/4
US630	2.	BAYONE	IN	IRR	3	7 3/4
US631	0	T	S	G	3	7 3/4
US632	0.	BAYONE	IN	IRR	4 1/2	8 1/2
US633	5	T	S	G	4 1/2	8 1/2
US634	0.	BAYONE	IN	IRR	4 1/2	8 1/2
US635	7	T	S	G	4 1/2	8 1/2
US636	1.	BAYONE	IN	IRR	4 1/2	8 1/2
US637	0	T	S	G	5	9
US638	1.	BAYONE	IN	IRR	5	9
US639	5	T	S	G	5	9
US640	2.	BAYONE	IN	IRR	5	9
US641	0	T	S	G	5	9
US532	0.	BAYONE	IN	IRR	6	10
US642	5	T	S	G	6	10
US643	0.	BAYONE	IN	IRR	6	10
US174	7	T	S	G	6	10
US533	1.	BAYONE	IN	IRR	6	10



US174-US533

CUSHING Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US369	0.	BAYONE	IN	IRR	3 1/2	7 1/2
US190	0.	BAYONE	IN	IRR	3 1/2	7 1/2
US191	0.	BAYONE	IN	IRR	3 1/2	7 1/2
US370	0.	BAYONE	IN	IRR	3 1/2	7 1/2



US190-US370



US190-US370



Electrosurgery

Irrigating Bipolar Forceps

HARDY Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US366	0.	BAYONE	INS	IRR	4 1/2	8 1/2
US382	5	T	INS	G	4 1/2	8 1/2
US385	1.	BAYONE	INS	IRR	4 1/2	8 1/2
US192	0	T	INS	G	4 1/2	8 1/2
US193	1.	BAYONE	INS	IRR	6	10
US194	5	T	INS	G	6	10
US195	2.	BAYONE	INS	IRR	6	10
US558	0	T	NON-INS	G	4 1/2	8 1/2

US192-US195, US382, US385, US558

US192-US195, US382, US385, US558

5 T G
 1. BAYONE IRR
 0. T 1/2 G
 1. BAYONE IRR
 5 T G
 0. BAYONE IRR
 5 T G

SCOVILLE Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US367	1.	BAYONE	INS	IRR	2 3/4	7 3/4
US556	5	T	NON-INS	G	2 3/4	7 3/4

US367, US556

1. BAYONE IRR
 5 T G
 1/2

GERALD Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US368	1.	BAYONE	INS	IRR	3 1/2	7 1/2
US557	0	T	NON-INS	G	3 1/2	7 1/2

US368, US557

1. BAYONE IRR
 0 T G
 1/2

MALIS Style Titanium Irrigating Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US360	STR	1.	BAYONE	IN	IRR	3 3/4	8 3/4
US361	STR	5	T	S	G	3 3/4	8 3/4
US362	STR	0.	BAYONE	IN	IRR	3 3/4	8 3/4
US363	ANG-UP	1.	T BAYONE	S IN	G IRR	3 3/4	8 3/4
US372	ANG-DOWN	0.5	T	S	G	3 3/4	8 3/4



US360-US362



US360-US362



US363, US372



US363, US372

MALIS Style Rose Gold Irrigating Bipolar Forceps

Catalog Number	Tip	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US166	STR	0.	BAYONE	IN	IRR	3 3/4	8 3/4
US167	STR	5	T	S	G	3 3/4	8 3/4
US168	STR	1.	BAYONE	IN	IRR	3 3/4	8 3/4
US169	ANG-UP	1.	T BAYONE	S IN	G IRR	3 3/4	8 3/4
US170	ANG-DOWN	0.5	T BAYONET BAYONE	S IN	G IRR	3 3/4	8 3/4



US166-US168



US169-US170





Electrosurgery

Irrigating Bipolar Forceps

Irrigating Bipolar Forceps, Stainless Steel

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
GK838	0.	BAYONE	IN	IRR	2	6 1/2
R	7	T	S	G	2 3/4	7
GK840	0.	BAYONE	IN	IRR	2 3/4	7
R	4	T	S	G	2 3/4	7
GK841	0.	BAYONE	IN	IRR	3 1/2	8
R	7	T	S	G	3 1/2	8
GK842	1.	BAYONE	IN	IRR	3 1/2	8
R	3	T	S	G	4 1/2	9
GK843	0.	BAYONE	IN	IRR	4 1/2	9
R	4	T	S	G	4 1/2	9
GK844	0.	BAYONE	IN	IRR		
R	7	T	S	G		
GK845	1.	BAYONE	IN	IRR		
R	3	T	S	G		
GK846	0.	BAYONE	IN	IRR		
R	4	T	S	G		
GK847	0.	BAYONE	IN	IRR		
R	7	T	S	G		
GK848	1.	BAYONE	IN	IRR		
R	3	T	S	G		



Electrosurgery

Disposable Bipolar Forceps

Stainless Steel Disposable Bipolar Forceps, Box of 10

Catalog Number	Tip Size (mm)	Shape	Style	Insulation	Irrigation	Working Length (in.)	Length (in.)
US153SU	1.5	BAYONET	SCOVILLE	INS	NON-IRRIG	2 3/4	7 3/4
US567SU	2.0	BAYONET	CUSHING	INS	NON-IRRIG	2 3/4	7 3/4
US562SU	0.7	BAYONET	CUSHING	INS	NON-IRRIG	3 1/2	7 1/2
US560SU	1.0	BAYONET	GERALD	INS	NON-IRRIG	3 1/2	7 1/2
US154SU	0.5	BAYONET	HARDY	INS	NON-IRRIG	4 1/2	8 1/2
US155SU	1.0	BAYONET	HARDY	INS	NON-IRRIG	4 1/2	8 1/2
US156SU	1.5	BAYONET	HARDY	INS	IRRIG	4 1/2	8 1/2
US157SU	2.0	BAYONET	HARDY	INS	IRRIG	4 1/2	8 1/2



US153SU-US567SU



US153SU-US567SU



US158SU



US158SU

ADSON Disposable Bipolar Forceps, Stainless Steel, Box of 10

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US158SU	1.0	STR	INS	NON-IRRIG	5

Disposable Jewelers Bipolar Forceps, Stainless Steel, Box of 10

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US568SU	0.4	STR	NON-INS	NON-IRRIG	4
US569SU	0.4	STR	INS	NON-IRRIG	4



US568SU-US569SU



ELECTRO SURGERY

INCHES

Electrosurgery

Disposable Bipolar Forceps



Rose Gold Disposable Jewelers Bipolar Forceps, Box of 5

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US408S	0.	ST	IN	NON-IRRG	4
U	4	R	S	NON-IRRG	5
US409S	1.	ST	IN		
U	0	R	S		



US408SU-US409SU

Rose Gold Disposable Bipolar Forceps, Box of 5

Catalog Number	Tip Size (mm)	Shape	Style	Insulation	Irrigation	Working Length (in.)	Length (in.)
US565SU	0.7	STR	SEMKIN	IN	NON-IRRG		5 1/2
US559SU	1.5	BAYONET	SCOVILLE	IN	NON-IRRG	2 3/4	7 3/4
US201SU	0.5	BAYONET	CUSHING	S	NON-IRRG	3 1/2	7 1/2
US561SU	0.7	BAYONET	CUSHING	IN	NON-IRRG	3 1/2	7 1/2
US202SU	1.0	BAYONET	CUSHING	S	NON-IRRG	3 1/2	7 1/2
US203SU	1.5	BAYONET	CUSHING	IN	NON-IRRG	3 1/2	7 1/2
US566SU	2.0	BAYONET	CUSHING	S	IRRG	3 1/2	7 1/2
US204SU	0.5	BAYONET	HARDY	IN	NON-IRRG	4 1/2	8 1/2
US563SU	0.7	BAYONET	HARDY	S	IRRG	4 1/2	8 1/2
US205SU	1.0	BAYONET	HARDY	IN	NON-IRRG	4 1/2	8 1/2
US206SU	1.5	BAYONET	HARDY	S	IRRG	4 1/2	8 1/2
US564SU	2.0	BAYONET	HARDY	IN	NON-IRRG	4 1/2	8 1/2



US201SU-US564SU



US565SU

IRRG
NON-IRRG
S
IN
S
IN
S
IN
S

Electrosurgery

Disposable Bipolar Forceps

Rose Gold Disposable Irrigating Bipolar Forceps, Box of 5

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US410S	0.	BAYONE	IN	IRR	4 1/2	8 1/2
U	5	T	S	G	4 1/2	8 1/2
US411S	1.	BAYONE	IN	IRR	4 1/2	8 1/2
U	0	T	S	G		
US412S	1.	BAYONE	IN	IRR		
U	5	T 1/2	S	G		

US410SU-US412SU

Stainless Steel Disposable Irrigating Bipolar Forceps, Box of 5

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Working Length (in.)	Length (in.)
US523S	0.	BAYONE	IN	IRR	4 1/2	8 1/2
U	5	T	S	G	4 1/2	8 1/2
US524S	1.	BAYONE	IN	IRR	4 1/2	8 1/2
U	0	T	S	G		
US525S	1.	BAYONE	IN	IRR		
U	5	T 1/2	S	G		

US523SU-US525SU

Specialty Disposable Non-Stick Bipolar Forceps

Catalog Number	Tip Size (mm)	Shape	Insulation	Irrigation	Length (in.)
US501S	0.	BAYONE	IN	NON-IRRIG	7
U	5	T	S	NON-IRRIG	7
US502S	1.	BAYONE	IN	NON-IRRIG	7
U	0	T	S	NON-IRRIG	8
US503S	1.	BAYONE	IN	NON-IRRIG	8
U	5	T	S	NON-IRRIG	8
US504S	0.	BAYONE	IN	NON-IRRIG	9
U	5	T	S	NON-IRRIG	9
US505S	1.	BAYONE	IN	NON-IRRIG	9
U	0	T	S		
US506S	1.	BAYONE	IN		
U	5	T	S		
US507S	0.	BAYONE	IN		
U	5	T	S		
US508S	1.	BAYONE	IN		
U	0	T	S		
US509S	1.	BAYONE	IN		
U	5	T	S		

US501SU-US509SU

1/1
US501SU-US509SU






ELECTRO SURGERY

INCHES

Electrosurgery

Bipolar Connecting Cords

Bipolar Cords for BipoJET®

Catalog Number	New Cord	Pin Distance	Compatible Generators
GN145		22 MM	EMC DOLLEX ERBE Int'l MARTIN
GN085		29 MM	AESCULAP VALLEYLAB ERBE
GN146		38 MM	CODMAN CMC II / III ERBE MARTIN



GK195



GN219



GN224

LANDOLT Bipolar Cord and Adapters

Catalog Number	Pin Distance (mm)	Description
GK195	28.6	REUSABLE CORD F/ INSULATED LANDOLT FORCEPS
GN219	32	ADAPTER F/ CONMED
GN224	38	ADAPTER F/ CODMAN

ELECTRO SURGERY

MILLIMETERS

Electrosurgery

Bipolar Connecting Cords



US349SP

Dual Bipolar Cord, Sterile Disposable

Catalog Number	Packaging	Length (ft)
US349SP	10/BOX	12



US351

Round Bipolar Cord, Reusable

Catalog Number	Material	Length (ft)
US351	SILICONE	12



US359

Dual Bipolar Cord, Reusable


Catalog Number	Length (ft)
US359	12



INCHES

Electrosurgery

Bipolar Connecting Cords



Combined Bipolar and Irrigation Cord

Catalog Number	Pin Type	Packaging	Compatible w/
US352	DUAL	BOX OF 10	-
US353	ROUND	BOX OF 10	-
US355	DUAL	BOX OF 10	MALIS 1000 IRRIGATION SYSTEM
US356	ROUND	BOX OF 10	MALIS 1000 IRRIGATION SYSTEM

ELECTRO SURGERY

